










Size and weight chart *VINCENTevolution4*



VINCENTevolution4 SIZES


					
	140 mm 5.51"	150 mm 5.90"	160 mm 6.29"	170 mm 6.69"	180 mm 7.08"
	XS	S	M	L	XL

WRIST OPTIONS WITH LAMINATION RING/PLATE						
VINCENTwrist quicksnap VINCENTring_coaxial_6 	Total length	175 mm 6.88 "	185 mm 7.28 "	195 mm 7.67 "	205 mm 8.07 "	215 mm 8.46 "
	Total weight	approx. 461 g	approx. 491 g	approx. 516 g	approx. 546 g	approx. 556 g
VINCENTwrist quicksnap + flexion VINCENTring_coaxial_6 	Total length	194 mm 7.63 "	197 mm 7.75 "	214 mm 8.42 "	217 mm 8.54 "	234 mm 9.21 "
	Total weight	approx. 506 g	approx. 516 g	approx. 561 g	approx. 571 g	approx. 601 g
VINCENTwrist short VINCENTtranscarpal_radial_4 	Total length	151 mm 5.94 "	154 mm 6.06 "	171 mm 6.73 "	174 mm 6.85 "	191 mm 7.51 "
	Total weight	approx. 412 g	approx. 422 g	approx. 467 g	approx. 477 g	approx. 507 g
VINCENTwrist short + flexion VINCENTtranscarpal_radial_4 	Total length	170 mm 6.69 "	173 mm 6.81 "	190 mm 7.48 "	193 mm 7.59 "	210 mm 8.26 "
	Total weight	approx. 457 g	approx. 467 g	approx. 512 g	approx. 522 g	approx. 552 g

Size and weight chart VINCENTyoung3+



VINCENTyoung3+



WRIST OPTIONS WITH LAMINATION RING/PLATE		
<p><i>VINCENTwrist quicksnap</i> <i>VINCENTring_coaxial_6</i></p>	Total length	166 mm 6.53 "
	Total weight	approx. 391 g
<p><i>VINCENTwrist quicksnap + flexion</i> <i>VINCENTring_coaxial_6</i></p>	Total length	185 mm 7.28 "
	Total weight	approx. 436 g
<p><i>VINCENTwrist short</i> <i>VINCENTtranscarpal_radial_4</i></p>	Total length	142 mm 5.59 "
	Total weight	approx. 347 g
<p><i>VINCENTwrist short + flexion</i> <i>VINCENTtranscarpal_radial_4</i></p>	Total length	161 mm 6.33 "
	Total weight	approx. 387 g