






Size and weight chart VINCENTyoung3+



VINCENTyoung3+



WRIST OPTIONS WITH LAMINATION RING/PLATE

<p>VINCENTwrist quicksnap <i>Lamination ring koaxial 6</i></p> 	<p>Total length</p>	<p>166 mm 6.53 "</p>
	<p>Total weight</p>	<p>approx. 391 g</p>
<p>VINCENTwrist quicksnap + flexion <i>Lamination ring koaxial 6</i></p> 	<p>Total length</p>	<p>185 mm 7.28 "</p>
	<p>Total weight</p>	<p>approx. 436 g</p>
<p>VINCENTwrist short <i>Transcarpal radial 4</i></p> 	<p>Total length</p>	<p>142 mm 5.59 "</p>
	<p>Total weight</p>	<p>approx. 347 g</p>
<p>VINCENTwrist short + flexion <i>Transcarpal radial 4</i></p> 	<p>Total length</p>	<p>161 mm 6.33 "</p>
	<p>Total weight</p>	<p>approx. 387 g</p>