

GF. glove factory

GF. color glove

TEXTILE GLOVES FOR HAND PROSTHESES



GF glove factory GmbH | Breite Straße 155 | 76135 Karlsruhe

Tel.: +49 721 161 784 62 | info@glovesfactory.de

www.glovesfactory.de



The textile gloves GF. color gloves are made of highly elastic knitted fabric with polyamide and Lycra® fibres. Therefore, its especially suitable for high agile hand prosthesis with single finger movement. These can be worn over the prosthetic hand as protection against contamination. The two-piece glove set consists of a underglove and an unicolor overglove. This product line of GF. color gloves are specially tailored for the dimensions of VINCENT hand prostheses. Furthermore, special designs for partial hand systems and prostheses of other manufacturers are also possible.

The properties

- Color black (other colors on request).
- Material: knitted fabric made of polyamid and Lycra® fibres.
- Slip-resistant palm in honeycomb design.
- Handschuhsystem besteht aus einem Unterhandschuh und einem Oberhandschuh.
- Protection against splash water and dirt due to water-repellent impregnation.
- The overgloves are available in five lengths and different circumference.

The advantages

- The high elasticity of the textile glove enables unrestricted movements of the prosthesis.
- Minimal wrinkling, as the textile is pre-tensioned to the prosthesis.
- Easy glove change.
- Washable.

Additional options

- Individual adaptation of the gloves for partial hand prostheses.

Storage life

- Depending on the degree of use, the gloves can be worn for 2-8 weeks.

